

Speaker Ready Room

Problems

Problem	The SRR is not using MongoDB at all, need to change that to a local db
Solution	Since SRR will be used without internet, the MongoDB must be local on the same machine. But during the development, we could use the database on 210.245.150.19 instead.
Fixed	
Problem	Need to convert the Session data to Topic instead
Solution	
Fixed	

Current Topic object, could be found in folder /data/topic

```
{
  "id": "114",
  "startTime": "2015-04-01T11:00:00Z",
  "endTime": "2015-04-01T12:30:00Z",
  "category": "invited",
  "subspecialty": "Glaucoma",
  "type": "symposium",
  "theme": "EGS Symposium: Cataract Surgery in Glaucoma",
  "venueId": "3",
  "venue": "Yuxiu Hall",
  "code": "GL-SS06",
  "title": "When is Clear-Lens Phacoemulsification Indicated in Primary Angle-Closure",
  "duration": "12",
  "sequence": "1",
  "abstractId": "0",
  "speakerLastName": "Husain",
  "speakerFirstName": "Rahat",
  "speakerBarcode": "59190",
  "speakerEmail": "rahat.husain@sneec.com.sg"
}
```

New Topic Object

```
{
  "code": "GL-SS06",
  "category": "Invited",
  "startTime": "11:00",
  "endTime": "12:30",
  "date": "2015-04-01",
  "subspecialty": "Glaucoma",
  "type": "Symposium",
  "theme": "EGS Symposium: Cataract Surgery in Glaucoma",
}
```

```
"venueId": "3",
  "venue": "Yuexiu Hall",
  "_id": "114",
  "id": "114",
  "abstractId": 0,
  "name": "When is Clear-Lens Phacoemulsification Indicated in Primary
Angle-Closure",
  "duration": 12,
  "sequence": 1,
  "speaker": {
    "id": "59190",
    "apaoId": "59190",
    "firstName": "Rahat",
    "gender": "",
    "middleName": "",
    "lastName": "Husain",
    "title": "Dr.",
    "biography": [],
    "subspecialties": [],
    "institution": {
      "name": "",
      "position": "",
      "department": ""
    },
    "displayName": "Rahat HUSAIN"
  }
}
```

Problem

From:
<https://wiki.questwork.com/> - Questwork's Wiki

Permanent link:
https://wiki.questwork.com/doku.php?id=development:speaker_ready_room:start&rev=1455364860

Last update: 2016/02/13 12:01

